TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) EXTERNAL FIXATION SYSTEM 3	EN	TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) EXTERNT FIXERINGSSYSTEM 89	SV
TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) SISTEMA DI FISSAZIONE ESTERNA 11	IT	TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) ΣΥΣΤΗΜΑ ΕΞΩΤΕΡΙΚΗΣ ΣΤΑΘΕΡΟΠΟΙΗΣΗΣ 97	EL
TRUELOK HEXAPODAL TL-HEX (TL-HEX) SYSTÈME DE FIXATION EXTERNE 20	FR	TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) 創外固定システム 107	JP
TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) EXTERNES FIXATIONSSYSTEM 29	DE	TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) 外固定系統 115	CN
SISTEMA TL-HEX TRUELOK HEXÁPODO (TL-HEX) SISTEMA DE FIJACIÓN EXTERNA 38	ES	TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) ZEVNÍ FIXAČNÍ SYSTÉM 122	CS
TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) EKSTERNT FIKSATIONSSYSTEM 47	DA	TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) SYSTEM STABILIZACJI ZEWNĘTRZNEJ 130	PL
TL-HEX TRUELOK HEXAPOD -JÄRJESTELMÄ (TL-HEX) ULKONINEN KIINNITYSJÄRJESTELMÄ 55	FI	SISTEM TL-HEX TRUELOK HEXAPOD (TL-HEX) SISTEM ZA ZUNANJO FIKSACIJO 139	SL
TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) EKSTERNE FIKSERINGSSYSTEM 63	NO	TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) 외고정 시스템 147	KR
TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) EXTERN FIXATIESYSTEEM 71	NL		
SISTEMA TL-HEX TRUELOK HEXAPOD (TL-HEX) SISTEMA DE FIXAÇÃO EXTERNA 80	PT		

ORTHOFIX Srl

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External fixation system

RX Only

Important information for the surgeon performing the operation

Read instructions before use

General information

The TL-HEX TRUELOK HEXAPOD SYSTEM (TL-HEX) consists of rings, foot plates and struts to be used in conjunction with the TrueLok external fixation system. The struts are made of assembly of stainless steel, aluminum and plastic components, whereas the circular rings are made of aluminum and stainless steel (locking screw, stud).

This external fixation system is a metal bone fixation device for stabilization and correction during the normal healing process. The purpose of the TL-HEX System is to provide a solution for deformity correction and keeping the bone segments aligned during the healing process; it is not designed to support weight bearing loads. The device may only be applied by a physician with adequate surgical and post-operative management knowledge of the system and its relative mechanical and software limitations.

The selection of the correct frame and implant size for each patient is fundamental for the success of the procedure. The metal screws and wires are subject to repeated stress during use and their resistance will be limited by the size and shape of the bones being treated. Unless care is taken in patient selection, proper placement of the rings, struts, wires and screws, and careful postoperative management to minimize the stresses on the frame, metal fatigue may occur, with consequent breakage or bending before the healing process is complete, which may result in further injury or the need to remove the device prematurely.

While using the TL-HEX system for deformity correction, the surgeon should keep the following in mind:

- Patient weight: an obese or overweight patient places loads on the device that could cause breakage or bending.
- Patient occupation or activity: the risk of bending or breakage of an internal or external bone fixation device during post-operative rehabilitation may increase if the patient carries out activities that involve lifting or heavy muscular activity, as these movements subject the device to forces that could cause it to break.
- Mental condition of the patient: the risk of breakage of a fixation device is greater in older patients, or patients with mental deficiency, alcoholics or drug addicts or patients who, for other reasons, may ignore the necessary restrictions and precautions to be observed while using the device.

- Patient sensitivity to certain material: if it is suspected that the patient is sensitive to the used material, suitable tests must be carried out before selecting or implanting the material.
- Debilitated patients: debilitated persons who have difficulty in using weight support devices can be at risk during post-operative rehabilitation.

The TL-HEX software can be accessed through the website www.tlhex.com for use by surgeons. The surgeons should have adequate training on the use of the software and instruct the patient to follow the prescriptions correctly.

The patient or his/her caretaker should be able to understand and make the adjustments as indicated in the prescription provided by the surgeon.

INDICATIONS

The TL-HEX System is intended for limb lengthening by metaphyseal or epiphyseal distractions, fixation of open and closed fractures, treatment of non-union or pseudoarthrosis of long bones and correction of bony or soft tissue defects or deformities.

Within this range, indications include:

- Post-traumatic joint contracture which has resulted in loss of range of motion
- Fractures and disease which generally may result in joint contractures or loss of range of motion and fractures requiring distraction
- Open and closed fracture fixation
- Pseudoarthrosis of long bones
- Limb lengthening by epiphyseal or metaphyseal distraction
- Correction of bony or soft tissue deformities
- Correction of bony or soft tissue defects
- Joint arthrodesis
- Infected fractures or non unions

CONTRAINDICATIONS

The TL-HEX System is not designed or sold for any use except as indicated. Use of the system is contraindicated in the following situation:

- Patients with mental or physiological conditions who are unwilling or incapable of following postoperative care instructions.

WARNINGS

- 1 Each implantable device is disposable. After use, the system must be discarded as it could be damaged and cause problems if reused.
- 2 The fixation device must never be subjected to bending, cutting or scratching as these could reduce the fixation frame's resistance to stress and increase the risk of bending or breakage.
- 3 Specialized instrumentation is required to apply the fixation devices.
- 4 Instruments subjected to prolonged use or excessive force are susceptible to fracture and should be examined prior to surgery.
- 5 The patient must be informed that a second minor surgery for the removal of the fixation system is required.
- 6 All equipment should be examined carefully prior to use in order to ensure it is working properly. If it is believed that a component or tool is defective, damaged or suspicious, DO NOT USE IT.

- 7 The locking screws on the rings and the strut locking bolt are retained and should never be removed for any reason. Rings and struts must be cleaned and sterilized in the assembled state with the locking screws in the ring and strut locking bolt **only** untightened.
- 8 If the strut stud does not fully seat in the ring mounting hole, back out the locking screw to the stop.
- 9 The fixation must be applied at such a distance from the skin as to allow post-surgical swelling and cleansing, bearing in mind that the stability of the fixation device depends on the distance between it and the bone. Additional equipment might be required for the application and removal of the fixation devices, such as wire-cutters, hammers and electric drills.
- 10 Do not use Magnetic Resonance imaging in segments where a fixator has been applied. Check the status of the screws and the fixation at regular intervals.
- 11 Do not overtighten the stud locking screw in the ring and the strut locking bolt as this causes stripping of their threads. Final locking of the set screw of the TL-HEX rings is performed with a torque wrench, which must be turned in a clockwise direction only. A click indicates that the correct torque has been applied. Any attempt to unlock the any screw using the torque wrench will damage its gearing. The torque wrench is pre-set at a specific value. This value should be checked at least every year or any time the instrument becomes damaged, by returning it to the local authorized representative.
- 12 If choosing two 5/8 Rings, the openings can be orientated in the same direction (i.e. both open anteriorly, posteriorly, medially). The proposed solution can anyway be not really applicable on the patient due to the position of the struts that can interfere with the soft tissues. The surgeon must check the real feasibility of the frame before applying it on the patient. This can be anyway overcome by adding a 3/8 ring. This component allows in fact to transform a 5/8 into a full ring, so, using this technique it's possible to apply two 5/8 rings with the opening on the same side (very useful in trauma cases).
- 13 When attaching ultrashort struts to the ring, do not insert both central studs into the same tab since this might cause interference between the two struts during assembly and correction manouvers. Avoid horizontal placement of the struts.
- 14 If choosing 300mm rings, use the 450mm or 550mm kirschner wires.
- 15 The patient/caregiver must be instructed regarding the use and maintenance of the fixation device and care of the wire/screw sites.
- 16 The patient/caregiver must be instructed by the physician to:
 - Make the adjustments or get help in making the adjustments as needed
 - Identify on the prescription when to return for a strut change and for follow-up visits
 - Report if adjustment schedule cannot be met
 - Report any adverse or unexpected effects (strut breakage or disengagement, components damage, clip dislodgement, lost prescription)
- 17 The surgeon must evaluate the integrity of the construct at follow-up visits.
- 18 The fracture or bone gap must be checked periodically during treatment, making any necessary adjustments to the fixation device. An excessive or persistent gap can delay consolidation.
- 19 In patients undergoing callus distraction, the regenerate bone must be checked regularly and monitored radiologically.
- 20 Removal of the device: the surgeon should decide when it is time to remove the fixation device. Before final removal, temporary removal of a section of the frame is advised to check the strength of the healed regenerate segment.
- 21 The rings tray (30110086) includes 5 removable inserts: 4 of which can hold 2 full rings, 2 5/8 rings and 2 3/8 rings of sizes greater than 120 mm, while the top insert can hold 2 full rings, 2 5/8 rings and 2 3/8 rings of

100 and 120 mm sizes. This configuration guarantees a weight within the ANSI AAMI ST77 Standards. Any combination of removable inserts different from the one in the tray (301 I0086) might exceed the maximum allowed weight.

POTENTIAL ADVERSE EFFECTS

1. Loosening, bending, rupture or fracture of the fixation device(s) or loss of fixation or migration which may result in nerve, soft tissue, or organ damage, including perforation through the skin or other hemorrhaging.
2. Loss of anatomical position with non-union or malunions with rotation or angulation.
3. Corrosion with localized tissue reaction or pain. Infection, local or systemic.
4. Pain, discomfort or abnormal sensations of the nervous system resulting from the presence of the device.
5. Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction.
6. Bone loss or reduced bone density due to a reduction in the tension applied to the bone.
7. Nerve or vascular damage following the insertion of wires or screws.
8. Deep or surface infections of the bone screw site, osteomyelitis, septic arthritis, including chronic drainage of the bone screw sites following removal of the device.
9. Oedema or possible compartmental syndrome.
10. Joint constriction, subluxation, dislocation or loss of motor movement.
11. Premature bone callus consolidation during distraction.
12. Possible tension affecting the soft tissues and/or the fixation during manipulation of the callus (e.g. corrections of deformities and/or elongation).
13. Lack of satisfactory bone regeneration, the development of mal-union or nonunion.
14. Fracture of regenerated bone, or at the bone screw holes, following removal of the device.
15. Loosening or breakage of the bone screws.
16. Bone damage due to erroneous bone screw selection.
17. Bone deformities or talipes equinus.
18. The persistence or recurrence of the initial condition subject to treatment.
19. New surgery to replace a component or the entire fixation frame.
20. Abnormal growth cartilage development in skeletally immature patients.
21. Foreign body reactions due to bone screws or components of the fixation frame.
22. Tissue necrosis secondary to the insertion of bone screws.
23. Pressure on the skin caused by external components when the free space is insufficient.
24. Limb dysmetria.
25. Excessive surgical bleeding.
26. Intrinsic risks associated with anesthesia.
27. Unmanageable pain.
28. Secondary bone sequestration due to rapid perforation of the cortex with accumulation of heat and bone necrosis.
29. Vascular disorders, including thrombophlebitis, pulmonary embolism, wound hematoma, avascular necrosis.

Warning: This device is not approved for fixing or attachment by means of screws to posterior elements (pedicles) of the cervical, thoracic or lumbar spinal column.

Important

Not all surgical cases are resolved with a positive result. Additional complications may arise at any time due to

improper use, for medical reasons or as a result of failure of the device, requiring surgery in order to remove or replace the external fixation device.

Presurgical and surgical procedures, including knowledge of the surgical techniques and the correct selection and positioning of the external fixation devices, are important aspects for the successful application of Orthofix external fixation devices by the surgeon.

Proper selection of the patient and their ability to comply with the physician's instructions and follow the prescribed treatment regimen will have a significant influence on the results. It is important to conduct patient screening and select the optimal therapy, while bearing in mind the requirements and/or limitations in terms of physical and/or mental activity. Should a candidate for surgery have any contraindications or be predisposed to any contraindications, DO NOT USE any Orthofix external fixation device.

TL-HEX Software

The Web-Based TL-HEX, TrueLok Hexapod System, Software aids the surgeon in his/her use of the TL-HEX product and it can be accessed at: www.tlhex.com.

Please contact customer support tlhexcustomer@orthofix.com to update account information.

The TL-HEX is a circular external fixator based on Ilizarov principles. The working segment of this frame consists of a hexapod made up of six variable length struts. The relative strut lengths determine the position of the rings in space. Because the rings are attached to bone segments, their position indirectly determines the position of the bone segments. The SW is able to calculate strut length adjustments for surgeon's review and approval. This software needs three sets of parameters to perform the calculation (see TL-HEX Software User Guide, available on the TL-HEX website).

With these sets of parameters, the computer will be able to calculate appropriate strut adjustments for surgeon's review and approval.

Identifiable information should not be entered into the data fields.

The Software allows healthcare professionals to delete their own records.

For further information, see TL-HEX Software User Guide and TL-HEX General Principles for the Frame Assembly.

Warnings

- 1 TL-HEX should be used only after careful study of the below reference guides and after careful determination of appropriate parameter values.
- 2 Refer to the TL-HEX reference guides: TL-HEX General Principles for the Frame Assembly, and TL-HEX Software User Guide for more information related to use of hardware and software. These documents contain a complete description of hardware and software and important information, such as indications for use, contraindications, warnings, and precautions.
- 3 Improper use of the software could result in erroneous or inadvertent calculations.
- 4 The software must be used only with TL-HEX. The software is not compatible with devices of other manufacturers.

- 5 It is important to check that the software diagrams correspond to the deformity that is seen on the patient's x-rays and/or clinically.
- 6 The surgeon must carefully review struts length adjustment calculations made by the SW for correctness and make sure print out is readable before giving it to the patient.
- 7 When using the software, click "Update Views" after any changes to the parameters.

Precautions

Accurate inputs are critical for accurate results. Verify and double check all input parameters. Intraoperative placement of the TL-HEX Fixator according to preoperative plans is imperative to achieve predetermined results. If intraoperative conditions require a change to frame placement (eccentricity) or size (parameters), new strut lengths will be calculated by entering the new inputs into the program. Small changes may affect accuracy of outcome. The End of Correction Screen is intended to identify if any struts go out of range during the treatment schedule. The Report should be referenced for final strut position (acute and gradual) and strut type.

Security Precautions

User is advised to clear the browser history (temporary Internet files, cookies, etc.) after logging out of the TL-HEX application.

Software Requirements

Recommended Browser and Display Settings: The use of the TL-HEX software requires interaction with the surgeon to input radiographic measurements and other parameters, followed by a visual review of the resulting TL-HEX assembly.

To ensure the best possible use of the software, the following are recommended:

- For use in Europe select the "language" setting in the Internet Explorer browser for your country of origin. This provides the correct date settings.
- Display Settings: Screen resolution of 1024 x 768 or higher.
- Supported Browsers: Microsoft Internet Explorer Version 8 or 9.
- Internet Connection: High-speed connectivity is recommended.

Materials

The TL-HEX struts are made of assembly of stainless steel, aluminum and plastic components, whereas the TL-HEX rings and foot plates are made of aluminum and stainless steel (locking screw, stud).

STERILE and NON STERILE products

Orthofix supplies certain external fixation devices as STERILE and others as NON STERILE. It is recommended the product label be checked in order to determine the sterility of each device.

Sterile

The devices or kits supplied as STERILE carry a label indicating said status.

The contents of the package are STERILE unless the package is opened or damaged. Do not use the device if the package has been opened or damaged.

Non Sterile

Unless indicated otherwise, the components of Orthofix external fixation devices are provided as NON STERILE. Orthofix recommends that all non sterile components be properly cleaned and sterilized, following the

specific cleaning and sterilization procedures recommended. The integrity and performance of the product is only guaranteed in the case the package is undamaged.

Cleaning and Maintenance (Please see instructions PQ ISP)

Prior to use, any NON STERILE product should be cleaned using a mixture of 70% medical grade alcohol and 30% distilled water. After cleaning, the device and/or system is duly rinsed in sterile distilled water and dried using a clean non-woven fabric cloth. Alkaline cleaning solutions should not be applied to the aluminum components during cleaning.

Prior to sterilization, it is recommended that all components be inspected, since any damage to the surface of the metal can reduce the fatigue strength and the robustness and lead to corrosion. Should the components be damaged in any way, they must be replaced with new ones.

It is recommended the fixation device be assembled in order to ensure that all the components are present.

Struts and rings should not be disassembled for cleaning. The locking screws in the ring and the locking bolt on the struts should be left untightened for cleaning and sterilization.

ALL DEVICES LABELED “SINGLE-USE ONLY” MUST NOT BE REUSED. ORTHOFIX IS ONLY RESPONSIBLE FOR THE SAFETY AND EFFICACY OF SINGLE USE COMPONENTS WITH THE FIRST PATIENT.

The institution or the user will be responsible for any subsequent use of said devices.

Sterilization

The recommended sterilization cycle is:

Method	Cycle	Temperature	Exposure time
Steam	Pre-vacuum	132°-135°C [270°- 275°F]	Minimum 10 minutes

Wrap the tray before sterilization with adequate material in order to prevent contamination after sterilization.

Do not or include additional systems or instruments in the sterilization tray.

Sterility cannot be guaranteed if the sterilization tray is overloaded.

Orthofix is only responsible for the safety and efficacy in the case of first patient's use of the external fixation devices.

The responsibility for any subsequent use will be entirely that of the institution or the attending physician.

PLEASE NOTE

In the case where the normal sensitivity of the limb is altered, abating the usual proprioceptive feedback, the fixation system may be subjected to greater loads than usual. Under these circumstances, the patient will have to be informed of the risk of excessive stress on the fixation system, and the physician should be especially attentive regarding problems associated with overloading which might lead to loosening, bending or breakage of the components. In such situations, it is recommended the rigidity of the fixation system be increased compared to what is usually necessary.

Implant removal

The external fixation devices are designed to stabilize the fracture site during the rehabilitation and/or limb lengthening process. Once the limb rehabilitation and/or lengthening is complete, the devices must be removed. Removal may also be recommended in other cases, for example if there are signs of corrosion, reaction or

localized pain; implant breakage, curvature, fracture or loosening; or bone loss.

Ordering information

Contact your local Orthofix sales representative.

TL-HEX TRUELOK HEXAPOD SYSTEM is a registered trademark of Orthofix Srl.

Note: the physician is fully responsible for the correct selection of the patient, proper training, experience in selecting and applying TL-HEX, as well as for the selection of the subsequent post-operative procedures. For complaints or additional information regarding the product and its use, contact Orthofix Srl.

RISKS DUE TO THE RE-USE OF “SINGLE USE” DEVICE

IMPLANTABLE DEVICE*

The “SINGLE USE” implantable device* of Orthofix is identified through symbol  reported on the product label.

After the removal from the patient, the implantable device* has to be discarded.

The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*): Implantable device

Any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

NON IMPLANTABLE DEVICE

The “SINGLE USE” non implantable device of Orthofix is identified through symbol  reported on the label or indicated in the “Instructions For Use” supplied with the products. The re-use of “SINGLE USE” non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

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All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

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